



# maxxeus™ DBM

## Allograft DBM Putty

- 100% Allograft
- Moldable
- Resistant to irrigation
- All lots tested for osteoinductivity
- All donors meet or exceed FDA and AATB standards
- BMP tested
- Demineralization process validated for viral inactivation



Description	Item Code
0.5 cc	2014-40
1.0 cc	2015-40
2.5 cc	2016-40
5.0 cc	2017-40
10.0 cc	2018-40

## Putty

DBM putty is packaged Sterile R in a plunger. The putty is moldable, resistant to irrigation and provides the maximum amount of DBM.

## Osteoinductive

Demineralized bone matrix (DBM) contains a variety of biologically active bone morphogenetic proteins (BMPs), which are among the osteoinductive components of bone. More than 20 BMPs have so far been identified in human allograft bone tissue. A validated in vitro osteoinductive assay using ELISA BMP 2 concentration measurement method is used to verify BMP is present in our final DBM product.



### Center for Tissue, Innovation and Research

2900 College Drive  
Kettering, Ohio 45420-2972  
1-855-MAXXEUS (629-9387)